# Instructions for use

## LP Rotate™ with EVO™

Product number:


VLXExxyz (foot module)

LRPOxy (male pyramid module)





LRFOxy (female pyramid module)





*Life Without Limitations®*

.....	3
Ⓔ EN Instructions for use .....	5
Ⓔ DE Gebrauchsanweisung .....	8
Ⓔ FR Notice d'utilisation .....	11
Ⓔ ES Instrucciones para el uso .....	14
Ⓔ IT Istruzioni per l'uso .....	17
Ⓔ NO Bruksanvisning .....	20
Ⓔ DA Brugsanvisning .....	23
Ⓔ SV Bruksanvisning .....	26
Ⓔ NL Gebruiksaanwijzing .....	29
Ⓔ PT Instruções de Utilização .....	32
Ⓔ JA 取扱説明書 .....	35
Ⓔ ZH 中文说明书 .....	38
Ⓔ KO 사용 설명서 .....	41

## LP ROTATE™ FOOT SYSTEM (Figure A)

- Low profile foot module with adjustable shock absorber and rotation with progressive stiffening
- Build height: with male pyramid adapter 145 - 153mm (5 11/16" - 6" ) with female pyramid adapter 153 – 165mm (6" - 6 1/2")
- Male pyramid adapter
- Female pyramid adapter
- Maintenance free, no lubrication required
- Maximum weight: 147kg (325lbs)
- Available in sizes 22-30

These recommendations are meant as starting points only. You may have to move up or down the chart based on amputee's comfort, impact level, or compression element.

Rod Color	Weight	Category
Yellow	45 – 68 Kg, 100 – 150 lbs	1 – 3
Blue	69 – 100 Kg, 151 – 220 lbs	4 – 6
Red	101 – 147 Kg, 221 – 324 lbs	7 – 8

## Bench Alignment (Figure B&C)

- Fit foot with selected cover.
- Adjust prosthesis to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/extension and abduction/adduction.
- Divide the foot cover into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram.

## Initial Fitting

When the user is standing, with equal weight on both limbs and full weight transferred onto the prosthesis, the LP Rotate™ should compress a maximum of 1-3mm. The compression resistance can be adjusted by either changing the inner rod or by turning the adjustment plug through the top of the male pyramid using 6mm hex key. It is sometimes necessary to upgrade the stiffness of the rod to control the vertical displacement and allow increased heel to toe function.

## Rod Adjustments (Figure D)

Compression resistance can be adjusted by either changing the compression element or by turning the adjustment plug (a) through the top of the male or female pyramid adapter using a 6mm hex key. Ensure adjustment plug (a) is in contact with inner rod. Turn the adjustment plug (a) CLOCKWISE for stiffer compression. COUNTERCLOCKWISE for a softer compression.

## Rod Replacement

The LP Rotate™ shock module contains an inner rod, which varies in density in order to control the vertical travel. This inner rod can be exchanged or replaced to customize the amount of the travel for patients of different weights. Compression resistance can be adjusted by following the rod adjustment.

1. Unscrew male or female pyramid adapter
2. Unscrew adjustment plug using a 6mm hex
3. Replace the inner rod
4. Ensure adjustment plug is in contact with inner rod.

5. Torque the pyramid adapter (**Figure E**) to specifications.
6. The pyramid clamp must be positioned anterior (never medial or lateral).

### **Dynamic Alignment**

The heel stores energy after initial contact, slowly releasing it at midstance. The forward momentum generated by this action results in the toe being loaded for optimum energy release at terminal stance. The heel to toe action can be influenced by:

- A-P positioning of foot
- Dorsi-Plantarflexion
- Heel Stiffness
- Adjusting the inner rod
- A soft rod can lead to the foot feeling too stiff. Upgrading the rod will transform the flexibility of the foot.
- Check for shock module action and assess vertical displacement.

### **HEEL STIFFNESS**

#### **I. Heel too soft**

##### **Symptoms**

- Foot will come to flat position too early (amputees will feel like sinking into a hole).
- Extra energy is required to climb up over the toe.
- Toe will feel too stiff.
- Knee may hyperextend.

##### **Action**

- Add polyurethane heel wedge.
- Shift socket anteriorly (or foot posterior).

#### **II. Heel too hard**

##### **Symptoms**

- Rapid heel to toe movement
- At heel strike the amputee has poor control of his prosthesis
- Minimal energy return feeling
- Knee may get unstable

##### **Action**

- Replace heel wedge with a smaller one
- Shift socket posteriorly (or foot anteriorly).

### **Heel Wedges CM150011**

The wedge can influence the heel to toe function. Start altering the heel resistance to improve foot response by adding a heel wedge. The 2°, 4° & 6° wedges can be interchanged to customize the stiffness and achieve the desired functional characteristics.

- 1 2°– yellow
- 2 4°– red
- 3 6°– blue

#### **I. Temporary Wedge Placement**

- Cut the wedge to the width of the foot module.
- Roughen the upper and lower surface of the wedge with abrasive paper.
- Place the wedge in the angle of heel and foot module.
- Secure in position bond to upper surface with tape wrapped around the foot module. (**Figure F**)

#### **II. Permanent Wedge Placement**

- Apply adhesive on the upper side of the wedge only.
- Locate in the foot/heel junction and position before adhesive sets.
- For split toe feet install the heel wedge then remove a thin slice in the middle by cutting with a sharp knife through the split in the carbon foot module. (**Figure G**)

Instant adhesive is necessary to bond the urethane heel wedge. The adhesive cures in 15-20 seconds. For removal the adhesive may be softened by soaking in acetone or cyanoacrylate adhesive remover.

### **Heel Divider**

The heel divider should be placed in the posterior of the splitted heel part. To secure its position a drop of instant adhesive can be used on one side.

### **Flex-Foot Sock with EVO**

The toe end of the sock is divided by a thread and fits within the toe split. Upon fitting the sock, hold the top of the toe thread and slide the sock into the toe split. The thread of the sock should rest on top of foot module and not underneath the toe split. The sock should be loose from the heel area and secured proximally with the cable tie provided.

### **Shoe Horn**

When removing and installing the cover use shoehorn to avoid damage to foot or cover. **(Figure H)**

**Note:** Tampering with the bolts that attach the shock unit with the foot module, will void the warranty of the product.

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**EN – Caution:** Össur products and components are designed and tested to ISO 10328. Compatibility and compliance with this standard is achieved only when Össur products and components are used with other recommended Össur or authorized components. If un-usual movement or product wear is detected in a structural part of a prosthesis at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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### **Össur Iberia S.A.**