Flex-Foot Axia®




Product number: LAPOXYYO




Life Without Limitations®

.....	3
(EN) Instructions for use.....	5
(DE) Gebrauchsanweisung.....	7
(FR) Notice d'utilisation.....	9
(ES) Instrucciones para el uso.....	11
(IT) Istruzioni per l'uso.....	13
(DA) Brugsanvisning.....	15
(SV) Bruksanvisning.....	17
(NL) Gebruiksaanwijzing.....	19
(PT) Instruções de Utilização.....	21
(JA) 取扱説明書.....	23

A**B****C**

Flex-Foot Axia® (Figure A)

Indicated for low to moderate impact levels only

- Maximum user weight 140 kg (308lbs)
- Heel height 10mm (3/8")
- Available in sizes 23-30
- Clearance 66 mm – 77 mm (2 5/8"-3")

Bench Alignment (Figures B+C)

- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/ extension and abduction/adduction.
- Divide the foot cover into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram (the guide inside of the foot cover matches the one third portion).

Dynamic Alignment

The tapered heel part with the compressed tarsal core stores energy at heel strike, then slowly releases it to assist tibial progression to midstance. The forward momentum generated by this action results in the toe being loaded for optimum energy release at terminal stance. The heel to toe action can be influenced by

- A-P positioning of foot
- Dorsi-Plantarflexion
- Heel stiffness adjustments

Heel Stiffness

Heel too soft

Symptoms

- Foot will come to flat position too early (amputee will feel like sinking into a hole)
- Extra energy is required to climb up over the toe
- Toe will feel too stiff.
- Knee may hyperextend.

Action

- Shift socket anterior (or foot posterior)
- Add inserts

Heel too hard

Symptoms

- Rapid heel to toe movement
- At heel strike amputee has poor control of his prosthesis
- Minimal energy return feeling
- Knee may become unstable

Action

- Shift socket posterior (or foot anterior)
- Remove the inserts

Heel Inserts (Figure D)

- A insert kit is available for the LAPOXYZZ.
- Instant adhesive is necessary to bond the heel inserts.
- The adhesive cures in 15-20 seconds.
- Care should be taken that it does not come into contact with skin.

Guidelines

- To insert: Put in the insert into the rib, furthest back on the tarsal core. Repeat on the other side (**Figure E**).

- To remove: Grab the exposed end of the insert with pliers and pull to remove.
- To attach permanently: Trim the insert to the width of the keel tarsal core.
- Apply a drop of Loctite 495 to the bottom of the rib and put in the insert.

Flex-Foot Sock FCX63006, FCX63007

To protect the cover and avoid noise a Flex-Foot sock is fitted over the foot module before installing the cover. The sock should be pressed between the toes, be loose in the heel area and secured proximally with the cable tie provided.

Shoe Horn CM100068 (Figure F)

When removing and installing the cover use shoehorn to avoid damage to foot or cover.

Note: Do not expose Tarsal Core to heat, abrasive materials or excessive moisture as it may cause damage.

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EN - Caution: Össur products and components are designed and tested to ISO 10328. Compatibility and compliance with this standard is achieved only when Össur products and components are used with other recommended Össur or authorized components. If un-usual movement or product wear is detected in a structural part of a prosthesis at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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Össur Iberia S.A.